+

+

+

D+WOUND SOLUTION

+

+

+

D+WOUND SOLUTION

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01. 당뇨병 DIABETIC FOOT ULCER | DIRECT CODING

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D+SOLUTION

- ▶ No debridement
- ▶ Moisturizing with Easydew
- ▶ Revascularization
- ▶ Cushioning with EasyFoam™

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D+SOLUTION

- ▶ Complete debridement
- ▶ CuraVAC® Silver
- ▶ basic fibroblast growth factor (bFGF)
- ▶ CGDerm™, CGCRYODERM®, CGDerm™ powder

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D+SOLUTION

- ▶ Complete debridement
- ▶ Systemic antibiotics
- ▶ CuraVAC® Silver
- ▶ bFGF
- ▶ CGDerm™, CGCRYODERM®, CGDerm™ powder

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D+SOLUTION

- ▶ EasyDERM™ Plus or EasyFoam™
- ▶ Easyef Ointment

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D+SOLUTION

- ▶ Confirm blood flow
- ▶ Complete debridement
- ▶ CuraVAC® Silver
- ▶ CGDerm™, CGCRYODERM®
- ▶ CGDerm™ powder
- ▶ bFGF, Easyef®, Easyef Ointment

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D+SOLUTION

- ▶ Evaluation of vascular status
- ▶ Complete debridement
- ▶ CuraVAC® Silver
- ▶ CGDerm™, CGCRYODERM®
- ▶ CGDerm™ powder
- ▶ bFGF, Easyef®, Easyef Ointment

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

D.1**Dry****I.0****No infection****R.1****Ischemic****E.0****No exudate****C.1****Chronic****T.0****Closed**

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OI. 당뇨병 DIABETIC FOOT ULCER | EXPERT OPINION

Flow Chart for DF according to **DIRECT**

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Topical agent of D+Soulution

Easyef®	Recombinant EGF
bFGF	Recombinant FGF
Easyef Ointment	the Ointment containing recombinant EGF
CGDerm™ paste	Acellular dermal matrix paste

OI. 당뇨병 DIABETIC FOOT ULCER | EXPERT OPINION

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E.2

E.3

How to convert vertical wound to horizontal by D+Solution

CuraVAC®	NPWT
bFGF	Recombinant FGF
CGDerm™	Acellular dermal matrix
CGDerm™ paste	Acellular dermal matrix paste

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01. 당뇨병 DIABETIC FOOT ULCER | REFERENCES

1. Jeffcoate WJ, Harding KG. Diabetic foot ulcers. Lancet 2003; 361(9368): 1545-51.
2. Wagner FW, Jr. The dysvascular foot: a system for diagnosis and treatment. Foot Ankle 1981; 2(2): 64-122.
3. Lavery LA, Armstrong DG, Harkless LB. Classification of diabetic foot wounds. J Foot Ankle Surg 1996; 35(6): 528-31.
4. Schaper NC. Diabetic foot ulcer classification system for research purposes: a progress report on criteria for including patients in research studies. Diabetes Metab Res Rev 2004; 20 Suppl 1: S90-5.
5. Hong JP, Jung HD, Kim YW. Recombinant human epidermal growth factor (EGF) to enhance healing for diabetic foot ulcers. Ann Plast Surg 2006; 56(4): 394-8.
6. Lee HB, Hong JP. Using recombinant human epidermal growth factor with standard of good care for diabetic foot ulcers in inoperable patients. J Korean Wound Management Soc 2008; 4(1): 16-9.

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D+SOLUTION

- ▶ Surgical or Autolytic or Enzymatic debridement
- ▶ CuraVAC® + Easyef® or bFGF

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D+SOLUTION

- ▶ Adequate Incision and Drainage and/or systemic antibiotics
- ▶ Surgical or Autolytic or Enzymatic debridement
- ▶ CuraVAC® + Easyef® or bFGF

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
PRESSURE ULCER ALGORITHM

1st

Adequate Incision and Drainage plus systemic antibiotics

2st

Surgical debridement

3rd

CuraVAC® + bFGF

4th

flap surgery

D.1

Dry

I.0

No-Infection

R.1

Ischemic

E.1

Light

C.1


Chronic

T.1

Open

1st

Surgical debridement


2st

Enzymatic debridement (EasyDERM™ Plus)

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1st

**Adequate incision and drainage plus
systemic antibiotics.**

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Stage 1	
Description	Skin intact but reddened for more than 1 hour after relief of pressure
D+WOUND SOLUTION	Easyef Ointment / EasyDERM™ Plus / EasyFoam™

Stage 2	
Description	Blister or other break in dermis and/or infection
D+WOUND SOLUTION	<p>1) No infection or infection control - EasyDERM™ Plus or EasyFoam™</p> <p>2) Infection</p> <ul style="list-style-type: none"> - Assess Clinical Sign and Symptoms of infection (purulent exudate and/or elevated temperature and/or peripheral induration or edema) - Adequate Incision and Drainage plus systemic antibiotics

Stage 3	
Description	Subcutaneous destruction into muscle and/or infection
D+WOUND SOLUTION	<p>Stage2 management Plus</p> <p>1) No infection or infection control</p> <ul style="list-style-type: none"> - EasyDERM™ Plus or EasyFoam™ - Surgical or Autolytic or Enzymatic debridement and CuraVAC® + Easyef® or bFGF <p>2) Wound Volume</p> <ul style="list-style-type: none"> - Small – EasyDERM™ Plus or EasyFoam™ - Large – Skin graft or Flap Surgery

Stage 4	
Description	Involvement of bone or joint and/or infection
D+WOUND SOLUTION	<p>Stage3 management Plus</p> <p>Surrounding Skin Condition</p> <p>1) Non-ischemic - EasyDERM™ Plus or EasyFoam™</p> <p>2) Ischemic</p> <ul style="list-style-type: none"> - Surgical or Autolytic or Enzymatic debridement & CuraVAC® + Easyef® or bFGF - Flap Surgery

O2. 욕창 | PRESSURE ULCER | REFERENCES

1. Reuler JB, Cooney TG. The pressure sore: pathophysiology and principles of management. Ann Intern Med 1981; 94(5): 661-6.
2. Shai A, Maibach HI. Wound healing and ulcers of the skin : diagnosis and therapy : the practical approach. Berlin ; New York: Springer; 2005.
3. Saha S, Smith MEB, Totten A, et al. Pressure Ulcer Treatment Strategies: Comparative Effectiveness. Rockville MD; 2013.
4. Livesley NJ, Chow AW. Infected pressure ulcers in elderly individuals. Clin Infect Dis 2002; 35(11): 1390-6.
5. de Laat EH, Scholte op Reimer WJ, van Achterberg T. Pressure ulcers: diagnostics and interventions aimed at wound-related complaints: a review of the literature. J Clin Nurs 2005; 14(4): 464-72.
6. Myers BA. Wound management : principles and practice. 2nd ed. Upper Saddle River, N.J.: Pearson/ Prentice Hall; 2008.

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
03. 외상 TRAUMATIC WOUND | REFERENCES

1. Leaper DJ. Traumatic and surgical wounds. BMJ 2006; 332(7540): 532-5.
2. Scherer SS, Pietramaggiori G, Mathews JC, Prsa MJ, Huang S, Orgill DP. The mechanism of action of the vacuum-assisted closure device. Plast Reconstr Surg 2008; 122(3): 786-97.
3. Scherer SS, Pietramaggiori G, Mathews JC, Orgill DP. Short periodic applications of the vacuum-assisted closure device cause an extended tissue response in the diabetic mouse model. Plast Reconstr Surg 2009; 124(5): 1458-65.
4. Franczyk M, Lohman RF, Agarwal JP, Rupani G, Drum M, Gottlieb LJ. The impact of topical lidocaine on pain level assessment during and after vacuum-assisted closure dressing changes: a double-blind, prospective, randomized study. Plast Reconstr Surg 2009; 124(3): 854-61.
5. Neligan P, Warren RJ, Van Beek A. Plastic surgery. 3rd ed. London ; New York: Elsevier Saunders; 2013.
6. Suissa D, Danino A, Nikolis A. Negative-pressure therapy versus standard wound care: a meta-analysis of randomized trials. Plast Reconstr Surg 2011; 128(5): 498-503.

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O4. 정맥성궤양 VENOUS ULCER | EXPERT OPINION

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O4. 정맥성궤양 | VENOUS ULCER | REFERENCES

1. O'Meara S, Martyn-St James M. Alginate dressings for venous leg ulcers. Cochrane Database Syst Rev 2013; 4: 1-74.
2. Meissner MH, Gloviczki P, Bergan J, et al. Primary chronic venous disorders. J Vasc Surg 2007; 46 Suppl S: 54S-67S.
3. Cronenwett JL, Johnston KW, Rutherford RB, Society for Vascular Surgery (U.S.). Rutherford's vascular surgery. 7th ed. Philadelphia, PA: Saunders/Elsevier; 2010.
4. Sabiston DC, Townsend CM. Sabiston textbook of surgery : the biological basis of modern surgical practice. 19th ed. Philadelphia, PA: Elsevier/Saunders; 2012.
5. Leopardi D, Hoggan BL, Fitridge RA, Woodruff PW, Maddern GJ. Systematic review of treatments for varicose veins. Ann Vasc Surg 2009; 23(2): 264-76.
6. Porter JM, Moneta GL. Reporting standards in venous disease: an update. International Consensus Committee on Chronic Venous Disease. J Vasc Surg 1995; 21(4): 635-45.
7. Eklof B, Rutherford RB, Bergan JJ, et al. Revision of the CEAP classification for chronic venous disorders: consensus statement. J Vasc Surg 2004; 40(6): 1248-52.
8. Bergan JJ, Schmid-Schonbein GW, Smith PD, Nicolaides AN, Boisseau MR, Eklof B. Chronic venous disease. N Engl J Med 2006; 355(5): 488-98.
9. Vasquez MA, Rabe E, McLafferty RB, et al. Revision of the venous clinical severity score: venous outcomes consensus statement: special communication of the American Venous Forum Ad Hoc Outcomes Working Group. J Vasc Surg 2010; 52(5): 1387-96.
10. CTD R, FAL D, GAF F. Hydrogel dressings for venous leg ulcers. The Cochrane Collaboration 2013: 1-15.
11. Hamdan A. Management of varicose veins and venous insufficiency. JAMA 2012; 308(24): 2612-21.
12. Robson MC, Cooper DM, Aslam R, et al. Guidelines for the treatment of venous ulcers. Wound Repair Regen 2006; 14(6): 649-62.
13. O'Meara S, Martyn-St James M. Foam dressings for venous leg ulcers. Cochrane Database Syst Rev 2013; 5: 1-114.
14. O'Donnell TF, Jr., Lau J. A systematic review of randomized controlled trials of wound dressings for chronic venous ulcer. J Vasc Surg 2006; 44(5): 1118-25.
15. Tang JC, Marston WA, Kirsner RS. Wound Healing Society (WHS) venous ulcer treatment guidelines: what's new in five years? Wound Repair Regen 2012; 20(5): 619-37.

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D+SOLUTION

- ▶ White-gray eschar
 - Chemical debridement
 - EasyDERM™ Plus for eschar softening
- ▶ Dark eschar
 - Surgical debridement
 - EasyDERM™ Plus for eschar softening

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D+SOLUTION

- ▶ 절개배농술(Incision and Drainage)
- ▶ Wound culture
- ▶ Topical/systemic antibiotics
- ▶ Debridement

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
D+SOLUTION

- ▶ Deep wound
→ CuraVAC® + CGDerm™
- ▶ Shallow wound → EasyFoam™

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D.1	I.0	R.1	E.0	C.1	T.0
Dry	No infection	Ischemic	No exudate	Chronic	Closed

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D.1**Dry****I.0****No infection****R.1****Ischemic****E.2****Intermediate****C.1****Chronic****T.1****Open**

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05. 동맥성궤양 ARTERIAL ULCER | EXPERT OPINION

1. Features of critical limb ischemia

① Physical Exam

- Dry skin, thickened nails, loss of hair, loss of subcutaneous fat or muscle atrophy
- Coolness to palpation
- Decreased or absent pulses
- Elevation pallor or dependent rubor
- Non-healing wound or ulcer, especially over bony prominences, distally, and on the plantar surface of the foot

② Non-invasive vascular

- Ankle-brachial index (ABI) ≤ 0.4
- Ankle systolic pressure ≤ 50 mmHg
- Toe systolic pressure ≤ 30 mmHg

Measures of skin microcirculation

(Capillary density ≤ 20 mm², Absent reactive hyperemia on capillary microscopy, TcPO₂ < 10 mmHg)

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- Barotraumas

05. 동맥성궤양 ARTERIAL ULCER | REFERENCES

1. A N, JR LJ. Critical limb ischemia: understanding the scope of the problem. Endovascular today. 2006;36-40.
2. Nelson EA, Bradley MD. Dressings and topical agents for arterial leg ulcers. The Cochrane database of systematic reviews 2007; (1); CD001836.
3. Fontaine R, Kim M, Kieny R. [Surgical treatment of peripheral circulation disorders]. Helv Chir Acta 1954; 21(5-6): 499-533.
4. Norgren L, Hiatt WR, Dormandy JA, et al. Inter-Society Consensus for the Management of Peripheral Arterial Disease (TASC II). Eur J Vasc Endovasc Surg 2007; 33 Suppl 1: S1-75.
5. Myers BA. Wound management : principles and practice. 2nd ed. Upper Saddle River, N.J.: Pearson/ Prentice Hall; 2008.
6. Shai A, Maibach HI. Wound healing and ulcers of the skin : diagnosis and therapy : the practical approach. Berlin ; New York: Springer; 2005.

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06. 화상 BURN | REFERENCES

1. Orgill DP, Ogawa R. Current methods of burn reconstruction. *Plast Reconstr Surg* 2013; 131(5): 827-36.
2. Cribbs RK, Luquette MH, Besner GE. Acceleration of partial-thickness burn wound healing with topical application of heparin-binding EGF-like growth factor (HB-EGF). *J Burn Care Rehabil* 1998; 19(2): 95-101.
3. Jahovic N, Guzel E, Arbak S, Yegen BC. The healing-promoting effect of saliva on skin burn is mediated by epidermal growth factor (EGF): role of the neutrophils. *Burns* 2004; 30(6): 531-8.
4. Hong JP, Kim YW, Jung HD, Jung KI. The effect of various concentrations of human recombinant epidermal growth factor on split-thickness skin wounds. *Int Wound J* 2006; 3(2): 123-30.
5. Alemdaroglu C, Degim Z, Celebi N, Sengezer M, Alomeroglu M, Nacar A. Investigation of epidermal growth factor containing liposome formulation effects on burn wound healing. *J Biomed Mater Res A* 2008; 85(1): 271-83.
6. Montesano R, Vassalli JD, Baird A, Guillemin R, Orci L. Basic fibroblast growth factor induces angiogenesis in vitro. *Proc Natl Acad Sci U S A* 1986; 83(19): 7297-301.
7. Haslik W, Kamolz LP, Manna F, Hladik M, Rath T, Frey M. Management of full-thickness skin defects in the hand and wrist region: first long-term experiences with the dermal matrix Matriderm. *J Plast Reconstr Aesthet Surg* 2010; 63(2): 360-4.
8. Shahrokhi S, Arno A, Jeschke MG. The use of dermal substitutes in burn surgery: acute phase. *Wound Repair Regen* 2014; 22(1): 14-22.
9. Akita S, Akino K, Imaizumi T, Hirano A. Basic fibroblast growth factor accelerates and improves second-degree burn wound healing. *Wound Repair Regen* 2008; 16(5): 635-41.

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O7. 수술봉합부위 SURGICAL SUTURE WOUND | DIRECT CODING

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
O7. 수술봉합부위 SURGICAL SUTURE WOUND | ALGORITHM

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O7. 수술봉합부위 SURGICAL SUTURE WOUND | EXPERT OPINION

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
O7. 수술봉합부위 SURGICAL SUTURE WOUND | REFERENCES

1. Kim EK, Hong JP. Efficacy of negative pressure therapy to enhance take of 1-stage allodermis and a split-thickness graft. *Ann Plast Surg* 2007; 58(5): 536-40.
2. Stannard JP, Robinson JT, Anderson ER, McGwin G, Jr., Volgas DA, Alonso JE. Negative pressure wound therapy to treat hematomas and surgical incisions following high-energy trauma. *J Trauma* 2006; 60(6): 1301-6.
3. Goldstein JA, Iorio ML, Brown B, Attinger CE. The use of negative pressure wound therapy for random local flaps at the ankle region. *J Foot Ankle Surg* 2010; 49(6): 513-6.
4. Suh HS, Lee AY, Park EJ, Hong JP. Negative pressure wound therapy on closed surgical wounds with dead space: animal study using a swine model. *Ann Plast Surg* 2014; doi: 10.1097/SAP.0000000000000231.
5. Scherer SS, Pietramaggiori G, Mathews JC, Prsa MJ, Huang S, Orgill DP. The mechanism of action of the vacuum-assisted closure device. *Plast Reconstr Surg* 2008; 122(3): 786-97.
6. Suissa D, Danino A, Nikolis A. Negative-pressure therapy versus standard wound care: a meta-analysis of randomized trials. *Plast Reconstr Surg* 2011; 128(5): 498e-503e.

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08. 열개창 WOUND DEHISCENCE | REFERENCES

1. Haynes AB, Weiser TG, Berry WR, et al. A surgical safety checklist to reduce morbidity and mortality in a global population. *N Engl J Med* 2009; 360(5): 491-9.
2. Horan TC, Andrus M, Dudeck MA. CDC/NHSN surveillance definition of health care-associated infection and criteria for specific types of infections in the acute care setting. *Am J Infect Control* 2008; 36(5): 309-32.
3. Clifford R. Wheelless J, M.D., Marcella L. Roenneburg MD. Atlas of Pelvic Surgery ; on-line edition. <http://www.atlasofpelvicsurgery.com/9AbdominalWall/5AbdominalWoundDehiscenceandEvisceration/cha9sec5.html>.
4. Mangram AJ, Horan TC, Pearson ML, Silver LC, Jarvis WR. Guideline for prevention of surgical site infection, 1999. Hospital Infection Control Practices Advisory Committee. *Infect Control Hosp Epidemiol* 1999; 20(4): 250-78.
5. Casey AL, Elliott TS. Progress in the prevention of surgical site infection. *Curr Opin Infect Dis* 2009; 22(4): 370-5.
6. Chong T, Sawyer R. Update on the Epidemiology and Prevention of Surgical Site Infections. *Curr Infect Dis Rep* 2002; 4(6): 484-90.
7. Fry DE. A systems approach to the prevention of surgical infections. *Surg Clin North Am* 2009; 89(2): 521-37.
8. McKibben L, Horan T, Tokars JL, et al. Guidance on public reporting of healthcare-associated infections: recommendations of the Healthcare Infection Control Practices Advisory Committee. *Am J Infect Control* 2005; 33(4): 217-26.
9. Services DoHaH. Action Plan to Prevent Healthcare-Associated Infections. <http://www.hhs.gov/ophs/initiatives/hai/infection.html>.

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Breast (62%) / Head and neck (24%) / Groin, genital and back (3%) / All other resions (8%).¹

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D+SOLUTION

- ▶ Autolytic or Enzymatic debridement

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09. 종양상처 CANCEROUS WOUND | REFERENCES

1. Sarabahi S, Tiwari VK, Bajaj SP. Principles and Practice of Wound Care. Jaypee Brothers Medical Pub; 2012.
2. Collier M. Management of patients with fungating wounds. Nurs Stand 2000;15:46-52
3. Froiland KG. Wound care of the advanced cancer patient. Hematol Oncol Clin North Am 2002;16:629-639.
4. Gerlach MA. Wound care issues in the patient with cancer. Nurs Clin North Am 2005;40:295-323.
5. Haubner F, Ohmann E, Pohl F, et al. Wound healing after radiation therapy: review of the literature. Radiat Oncol 2012;7:162.
6. O'Brien C. Malignant wounds: managing odour. Can Fam Physician 2012;58:272-274.

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
Taping

Laser Therapy

**Adjunctives :
Moisturizing
cream (Easydew)
/
Sun block**

**Skin Resurfacing
/Soft Tissue Filler
(Facetem™)**

**Scar Revision
/Scar Release
and Coverage**

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(2) Laser therapy (Secondary management)

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IO. 흉터 SCAR | REFERENCES

1. Mustoe TA, Cooter RD, Gold MH, et al. International clinical recommendations on scar management. *Plast Reconstr Surg* 2002; 110(2): 560-71.
2. Steinstraesser L, Flak E, Witte B, et al. Pressure garment therapy alone and in combination with silicone for the prevention of hypertrophic scarring: randomized controlled trial with intraindividual comparison. *Plast Reconstr Surg* 2011; 128(4): 306e-13e.
3. Stavrou D, Weissman O, Winkler E, et al. Silicone-based scar therapy: a review of the literature. *Aesthetic Plast Surg* 2010; 34(5): 646-51.
4. Mustoe TA. Evolution of silicone therapy and mechanism of action in scar management. *Aesthetic Plast Surg* 2008; 32(1): 82-92.
5. Akaishi S, Akimoto M, Hyakusoku H, Ogawa R. The tensile reduction effects of silicone gel sheeting. *Plast Reconstr Surg* 2010; 126(2): 109e-11e.

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
Product Publication List

1. Cho JH, Lew DH, Kim YS, Tark KC. Treatment of acute and chronic wound in immunodeficiency patients using Easyef® (EGF). *J Korean Soc Transplant* 2006; 20(1): 99–103.
2. Hong JP, Jung HD, Kim YW. Recombinant human epidermal growth factor (EGF) to enhance healing for diabetic foot ulcers. *Ann Plast Surg* 2006; 56(4): 394–8.
3. Jung K, Kim S, Moon SY, et al. Effects of recombinant human epidermal growth factor (rhEGF) on experimental radiation-induced oral mucositis in rats. *J Korean Soc Ther Radiol Oncol* 2006; 24(1): 67–76.
4. Kim HS, Kang KM, Lee Sw, Na JB, Chai GY. Effects of recombinant human epidermal growth factor on the proliferation and radiation survival of human fibroblast cell lines in vitro. *J Korean Soc Ther Radiol Oncol* 2006; 24(3): 179–84.
5. Jung K, Park HJ, Kim YH, Hong JP. Expression of tumor necrosis factor- α and cyclooxygenase-2 mRNA in porcine split-thickness wounds treated with epidermal growth factor by quantitative real-time PCR. *Res Vet Sci* 2007; 82(3): 344–8.
6. Lee SW, Jung KI, Kim YW, Jung HD, Kim HS, Hong JP. Effect of epidermal growth factor against radiotherapy-induced oral mucositis in rats. *Int J Radiat Oncol Biol Phys* 2007; 67(4): 1172–8.
7. Lee SW, Moon SY, Kim YH, Hong JP. The use of recombinant human epidermal growth factor to promote healing for chronic radiation ulcer. *Int Wound J* 2007; 4(3): 216–20.
8. Na JB, Kim HJ, Chai GY, et al. The effect of recombinant human epidermal growth factor on cisplatin and radiotherapy induced oral mucositis in mice. *J Korean Soc Ther Radiol Oncol* 2007; 25(4): 242–8.
9. Hong JP, Kim YW, Lee SK, Kim SH, Min KH. The effect of continuous release of recombinant human epidermal growth factor (rh-EGF) in chitosan film on full thickness excisional porcine wounds. *Ann Plast Surg* 2008; 61(4): 457–62.
10. Lee HB, Hong JP. Using recombinant human epidermal growth factor with standard of good care for diabetic foot ulcers in inoperable patients. *J Korean Wound Management Soc* 2008; 4(1): 16–9.
11. Lee KK, Jo HJ, Hong JP, et al. Recombinant human epidermal growth factor accelerates recovery of mouse small intestinal mucosa after radiation damage. *Int J Radiat Oncol Biol Phys* 2008; 71(4): 1230–5.
12. Yim H, Cho YS, Kim D-h, et al. The clinical effectiveness of Easyef® (Epidermal Growth Factor Spray) in pediatric burns. *J Korean Burn Soc* 2008; 11(2): 92–5.
13. Hong JP, Lee SW, Song SY, et al. Recombinant human epidermal growth factor treatment of radiation-induced severe oral mucositis in patients with head and neck malignancies. *Eur J Cancer Care* 2009; 18(6): 636–41.
14. Song JY, Lee SW, Hong JP, Chang SE, Choe H, Choi J. Epidermal growth factor competes with EGF receptor inhibitors to induce cell death in EGFR-overexpressing tumor cells. *Cancer Lett* 2009; 283(2): 135–42.
15. Tuyet HL, Nguyen Quynh TT, Vo Hoang Minh H, et al. The efficacy and safety of epidermal growth factor in treatment of diabetic foot ulcers: the preliminary results. *Int Wound J* 2009; 6(2): 159–66.
16. Wu HG, Song SY, Kim YS, et al. Therapeutic effect of recombinant human epidermal growth factor (RhEGF) on mucositis in patients undergoing radiotherapy, with or without chemotherapy, for head and neck cancer: a double-blind placebo-controlled prospective phase 2 multi-institutional clinical trial. *Cancer* 2009; 115(16): 3699–708.

Product Publication List

17. Cho KH, Kim YJ. Therapeutic Effect of the Recombinant Human Epidermal Growth Factor (rhEGF) in Pressure Ulcer. *J Korean Acad Rehabil Med* 2010; 34(3): 253–8.
18. Choi J, Moon SY, Hong JP, Song JY, Oh KT, Lee SW. Epidermal growth factor induces cell death in the absence of overexpressed epidermal growth factor receptor and ErbB2 in various human cancer cell lines. *Cancer Invest* 2010; 28(5): 505–14.
19. Kim JK, Kim CS, Ahn HJ, et al. Early recombinant human epidermal growth factor treatment recovers the irradiation-induced decrease of Na⁺ absorption prior to the definite histological mucositis. *Biomed Pharmacother* 2010; 64(9): 594–9.
20. Ryu SH, Kang KM, Moon SY, et al. Therapeutic effects of recombinant human epidermal growth factor (rhEGF) in a murine model of concurrent chemo- and radiotherapy-induced oral mucositis. *J Radiat Res* 2010; 51(5): 595–601.
21. Ryu SH, Kim YH, Lee SW, Hong JP. The preventive effect of recombinant human growth factor (rhEGF) on the recurrence of radiodermatitis. *J Radiat Res* 2010; 51(5): 511–7.
22. Hong JP. The use of recombinant human epidermal growth factor(rh-EGF, Nepidermin) to minimize scar formation. *J Wound Technol* 2010; 10: 22–5.
23. Kim YS, Lew DH, Tark KC, Rah DK, Hong JP. Effect of recombinant human growth factor against cutaneous scar formation in murine full-thickness wound healing. *J Korea Med Sci* 2010; 25(4): 589–96.
24. Kim EK, Li G, Lee TJ, Hong JP. The effect of human adipose-derived stem cells on healing of ischemic wounds in a diabetic nude mouse model. *Plast Reconstr Surg* 2011; 128(2): 387–94.
25. Lee JI, Kim KS, Oh BC, et al. Acute necrotic stomatitis (noma) associated with methicillin-resistant *Staphylococcus aureus* infection in a newly acquired rhesus macaque (*Macaca mulatta*). *J Med Primatol* 2011; 40(3): 188–93.
26. Hong JP, Park SW. The combined effect of recombinant human epidermal growth factor and erythropoietin on full-thickness wound healing in diabetic rat model. *Int Wound J* 2012; doi: 10.1111/j.1742-481X.2012.01100.x.
27. Ben-Nakhi M, Park EJ, Hong JP. The combined use of recombinant human epidermal growth factor (rh-EGF, Nepidermin) and fractional carbon dioxide laser in reducing cutaneous scar formation. *J Wound Technol* 2012; 15: 18–20.
28. Hong JP. Topical rh-EGF (easyef) for various wounds and scars. *J Wound Technol* 2013; 22.
29. Lee JH, Lee, SW, Hong JP, Shon MW, Ryu SH, Ahn SD. Foam dressing with epidermal growth factor for severe radiation dermatitis in head and neck cancer patients. *Int Wound J* 2014; doi: 10.1111/iwj.12317.

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Product Publication List

1. Kim EK, Hong JP. Efficacy of negative therapy to enhance take of 1-stage allodermis and a split-thickness graft. *Ann Plast Surg* 2007; 58: 536–40.
2. Lee KN, Seo DM, Hong JP. The effect and safety after extended use of continuous negative pressure of 75 mmHg over mesh and allodermis graft on open sternal wound from oversized heart transplant in a 3-month-old infant. *Int Wound J* 2010; 7(5): 379–84.
3. Lee SP, Ji SY, Lee JW, Yang WS. Efficacy of curavac system to enhance skin graft survival in chronic intractable wound. *J Korean Wound Management Soc* 2010; 6(1): 18–24.
4. Kim BJ, Suk JH, Jo AR, et al. Negative-pressure wound therapy using modified vacuum-assisted closure in patients with diabetic foot ulcers. *J Korean Diabetes* 2011; 12(2): 122–7.
5. Han WS, Kim K. Acute Postpneumectomy empyema with bronchopleural fistula treated with vacuum-assisted closure device. *Korean J Thorac Cardiovasc Surg* 2012; 45(4): 260–2.
6. Kim YH, Kim JS, Young KW, Lee HS, Hur JU, Lee KT. Treatment of a chronic gouty ulcer on the lateral malleolus with vacuum-assisted closure (VAC) therapy. *J Rheum Dis* 2012; 19(4): 230–3.
7. Lee KN, Ben-Nakhi M, Park EJ, Hong JP. Cyclic negative pressure wound therapy: an alternative mode to intermittent system. *Int Wound J* 2013; doi: 10.1111/iwj.12201.
8. Oh BH, Lee SH, Nam KA, Lee HB, Chung KY. Comparison of negative pressure wound therapy and secondary intention healing after excision of acral lentiginous melanoma on the foot. *Br J Dermatol* 2013; 168(2): 333–8.
9. Suh H, Hong JP. One stage allogenic acellular dermal matrices (ADM) and split-thickness skin graft with negative pressure wound therapy. *InTech* 2013; <http://dx.doi.org/10.5772/53304>.
10. Lee KT, Pyon JK, Lim SY, Mun GH, Oh KS, Bang SI. Negative-pressure wound dressings to secure split-thickness skin grafts in the perineum. *Int Wound J* 2014; 11(2): 223–7.
11. Suh H, Lee AY, Park EJ, Hong JP. Negative pressure wound therapy on closed surgical wounds with dead space: animal study using a swine model. *Ann Plast Surg* 2014; doi: 10.1097/SAP.0000000000000231.

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Product Publication List

1. Lee KN, Ben-Nakhi M, Park EJ, Hong JP. Cyclic negative pressure wound therapy: an alternative mode to intermittent system. *Int Wound J* 2013; doi: 10.1111/iwj.12201.
2. Suh H, Hong JP. One stage allogenic acellular dermal matrices (ADM) and split-thickness skin graft with negative pressure wound therapy. *InTech*. 2013; <http://dx.doi.org/10.5772/53304>.
3. Lee KT, Pyon JK, Lim SY, Mun GH, Oh KS, Bang SI. Negative-pressure wound dressings to secure split-thickness skin grafts in the perineum. *Int Wound J* 2014; 11(2): 223-7.

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Product Publication List

1. Kim EK, Hong JP. Efficacy of negative pressure therapy to enhance take of 1-stage allodermis and a split-thickness graft. Ann Plast Surg 2007; 58(5): 536-40.
2. Suh H, Hong JP. One stage allogenic acellular dermal matrices (ADM) and split-thickness skin graft with negative pressure wound therapy. InTech 2013; <http://dx.doi.org/10.5772/53304>.
3. Lee JH, Park KR, Kim TG, et al. A comparative study of CGCryoDerm and Alloderm in direct-to-implant immediate breast reconstruction. Arch Plast Surg 2013; 40(4): 374-79.

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Product Publication List

1. Hong JP, Kim YW, Jung HD, Jung KI. The effect of various concentrations of human recombinant epidermal growth factor on split-thickness skin wounds. Int Wound J 2006; 3(2): 123-30.
2. Hong JP. The use of recombinant human epidermal growth factor(rh-EGF, Nepidermin) to minimize scar formation. J Wound Technol 2010; 10: 22-5.
3. Kim YS, Lew DH, Tark KC, Rah DK, Hong JP. Effect of recombinant human epidermal growth factor against cutaneous scar formation in murine full-thickness wound healing. J Korean Med Sci 2010; 25(4): 589-96.

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Product Publication List

1. Ryu SH, Kim YH, Lee SW, Hong JP. The preventive effect of recombinant human growth factor (rhEGF) on the recurrence of radiodermatitis. *J Radiat Res* 2010; 51(5): 511–7.
2. Rhee DY, Park GH, Rho NK, Chang SE. The prediction and prevention of laser-induced post-inflammatory hyperpigmentation. *Korean J Cosmet Dermatol* 2012; 9:39–42.
3. Ben-Nakhi M, Park EJ, Hong JP. The combined use of recombinant human epidermal growth factor (rh-EGF, Nepidermin) and fractional carbon dioxide laser in reducing cutaneous scar formation. *J Wound Technol* 2012; 15:18–20.
4. Kong M, Hong SE. Topical use of recombinant human epidermal growth factor (EGF)-based cream to prevent radiation dermatitis in breast cancer patients: a single-blind randomized preliminary study. *Asian Pac J Cancer Prev* 2013; 14(8): 4859–64.
5. Yun WJ, Bang SH, Min KH, Kim SW, Lee MW, Chang SE. Epidermal growth factor and epidermal growth factor signaling attenuate laser-induced melanogenesis. *Dermatol Surg* 2013; 39(12): 1903–11.
6. Wu et al. The safety and efficacy of EGF-based cream for the treatment of radiotherapy-induced skin injury: Results from a multicenter observational study. (2014, submitted).
7. Lee JH, Lee SW, Hong JP, Shon MW, Ryu SH, Ahn SD. Foam dressing with epidermal growth factor for severe radiation dermatitis in head and neck cancer patients. *Int Wound J* 2014; doi: 10.1111/iwj.12317.

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APPENDICES

APPENDIX I

Classifications of Diabetic Foot Ulcer

1. Wagner classification

Grade	Lesion
0	No open lesions; may have deformity or cellulitis
1	Superficial diabetic ulcer (partial or full thickness)
2	Ulcer extension to ligament, tendon, joint capsule, or deep fascia without abscess or osteomyelitis
3	Deep ulcer with abscess, osteomyelitis, or joint sepsis
4	Gangrene localized to portion of forefoot or heel
5	Extensive gangrenous involvement of the entire foot

Wagner FW, Jr. The dysvascular foot: a system for diagnosis and treatment. Foot Ankle 1981; 2(2): 64-122.

2. UT (University of Texas) grading

		Grade			
		0	1	2	3
Grade	A	Pre or postulcerative lesion completely epithelialized	Superficial wound, not involving tendon, capsule or bone	Wound penetrating to tendon or capsule	Wound penetrating to bone or joint
	B	With infection	With infection	With infection	With infection
	C	With ischemia	With ischemia	With ischemia	With ischemia
	D	With infection and ischemia	With infection and ischemia	With infection and ischemia	With infection and ischemia

Lavery LA, Armstrong DG, Harkless LB. Classification of diabetic foot wounds. J Foot Ankle Surg 1996; 35(6): 528-31.

APPENDIX I

Classifications of Diabetic Foot Ulcer

3. PEDIS grading

Grade	Infection Severity	Clinical Manifestations	Treatment Parameters	Medications
1	Uninfected	Wound without purulence or inflammation	Outpatient	Topical antibiotics
2	Mild	≥2: purulence or erythema, pain, tenderness, warmth, or induration; cellulitis ≤2cm around ulcer; infection limited to skin/subcutaneous tissue; no other complications	Most not limb-threatening; most outpatient treatment	Cephalexin, trimethoprim-sulfamethoxazole(TMP-SMX), levofloxacin, amoxicillin-clavulanate, clindamycin
3	Moderate	Infection as above plus >1: cellulitis>2cm, streaking, deep tissue abscess, gangrene and with some life-threatening; involvement of muscle, tendon, joint, or bone	Most limb-threatening with some life-threatening; requires hospital treatment	TMP-SMX, amoxicillin-clavulanate, levofloxacin, ceftriaxone, linezolid, ertapenem, ticarcillin-clavulanate
4	Severe	Infection plus systemic toxicity or metabolic instability; fever, chills, tachycardia, hypotension, confusion, vomiting, severe hyperglycemia, acidosis, or azotemia	Life-threatening; requires hospital treatment	Imipenem-cilastatin, vancomycin-ceftazidime, levofloxacin-clindamycin, piperacillin-tazobactam, ticarcillin-clavulanate

Schaper NC. Diabetic foot ulcer classification system for research purposes: a progress report on criteria for including patients in research studies. *Diabetes Metab Res Rev* 2004; 20 Suppl 1: S90-5.

TIP

The Wagner classification is one of the most popular classification systems. The classification is simple and easy to use but only based on wound depth and necrosis, and thus it is limited to plan appropriate treatment strategies. The PEDIS grading is mainly used in the division of infectious disease. The UT classification is considered as the most appropriate system for wound specialists because it assesses ulcer depth and the presence of wound infection along with the presence of clinical signs of lower-extremity ischemia.

APPENDIX II

Classification of Pressure Ulcer

NPUAP (National Pressure Ulcer Advisory Panel) Classification

Stage. 1

Intact skin with non-blanchable redness of a localized area usually over a bony prominence. Darkly pigmented skin may not have visible blanching; its color may differ from the surrounding area.

Stage. 2

Partial thickness loss of dermis presenting as a shallow open ulcer with a red pink wound bed, without slough. May also present as an intact or open/ruptured serum-filled blister.

Stage. 3

Full thickness tissue loss. Subcutaneous fat may be visible but bone, tendon or muscle are not exposed. Slough may be present but does not obscure the depth of tissue loss. May include undermining and tunneling.

Stage. 4

Full thickness tissue loss with exposed bone, tendon or muscle. Slough or eschar may be present on some parts of the wound bed. Often includes undermining and tunneling.

APPENDIX III

Classification of Chronic Venous Disorders

CEAP classification

CEAP	Description
1. Clinical classification	
C0	no visible or palpable signs of venous disease
C1	telangiectasies or reticular veins
C2	varicose veins
C3	edema
C4a	pigmentation or eczema
C4b	lipodermatosclerosis or atrophie blanche
C5	healed venous ulcer
C6	active venous ulcer
2. Etiologic classification	
Ec	congenital
Ep	primary
Es	secondary (post-thrombotic)
En	no venous cause identified
3. Anatomic classification	
As	superficial veins
Ap	perforator veins
Ad	deep veins
An	no venous location identified
4. Pathophysiologic classification	
Pr	reflux
Po	obstruction
Pr.o	reflux and obstruction
Pn	no venous pathophysiology identifiable

Eklof B, Rutherford RB, Bergan JJ, et al. Revision of the CEAP classification for chronic venous disorders: consensus statement. J Vasc Surg 2004; 40(6): 1248-52.

APPENDIX IV

Classification of Peripheral Arterial Disease

Classification of Fontaine¹ and Rutherford²

Fontaine		Rutherford		
Stage	Clinical	Grade	Category	Clinical
I	Asymptomatic	0	0	Asymptomatic
IIa	Mild claudication	I	1	Mild claudication
IIb	Moderate-severe claudication	I	2	Moderate claudication
		I	3	Severe claudication
III	Ischemic rest pain	II	4	Ischemic rest pain
IV	Ulceration or gangrene	III	5	Minor tissue loss
		IV	6	Ulceration or gangrene

1. Fontaine R, Kim M, Kiény R. [Surgical treatment of peripheral circulation disorders]. *Helv Chir Acta* 1954; 21(5-6): 499-533.

2. Norgren L, Hiatt WR, Dormandy JA, et al. Inter-Society Consensus for the Management of Peripheral Arterial Disease (TASC II). *Eur J Vasc Endovasc Surg* 2007; 33 Suppl 1: S1-75.

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© 시지바이오, 2014